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K=A-B

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Real Time PCR Report

1 Laboratory equipment and reagents

1.1 Laboratory equipment and consumable

Equipment	Manufacturers	Model
High-speed-microtherm Homogenizer	Wanwu	KZ-III-F
Centrifuge	DragonLab	D3024R
Real time PCR System	ABI	Stepone plus
Clean bench	Suzhou Antai Air Tech Co., Ltd	SW-CJ-1FD
Ultramicro spectrophotometer	Thermo	NanoDrop2000
Standard Reagent Type Lab Ultra pure Water Purifier	Qingdao Fulum Technology Co., Ltd	FBZ2001-up-p
Centrifuge tube	Wanwu	
TIP	Wanwu	

1.2 Reagents

Reagents	Manufacturers	Order
RNA Extraction	Wanwu	G3013
Trichloromethane	Sinopharm Group Chemical Reagent Co., Ltd.	10006818
iso-Propyl alcohol	Sinopharm Group Chemical Reagent Co., Ltd.	80109218

Ethanol	Sinopharm Group Chemical Reagent Co., Ltd.	10009218
HyPure™Molecular Biology Grade Water	HyClone	SH30538.02
Wanwu®RT First Strand cDNA Synthesis Kit	Wanwu	G3330
2×SYBR Green qPCR Master Mix (High ROX) Primer	Wanwu Wanwu	G3322

2 Protocol:

2.1 Isolate RNA

2.1.1 Homogenization

2.1.1.1 Tissues: Homogenize tissue samples in 1 ml of RNA Extraction per 100 mg of tissue using power homogenizer. The sample volume should not exceed 10% of the volume of RNA Extraction used for the homogenization.

2.1.1.2 Cells grown in Monolayer: Rinse cell monolayer with ice cold PBS once. Lyse cells directly in a culture dish by adding 1 ml of RNA Extraction per 3.5 cm diameter dish. Pass the cell lysate several times through a pipette. Vortex thoroughly. The amount of RNA Extraction added is based on the area of the culture dish (1 ml per 10 cm²) and not on the number of cells present. An insufficient amount of RNA Extraction may result in DNA contamination of the isolated RNA.

2.1.1.3 Cells Grown in suspension: Spin cells for 10 min at 3000rpm. Remove media and resuspend cells in ice cold PBS. Pellet cells by spinning at 3000rpm for 10 min. Lyse cells with RNA Extraction by repetitive pipetting. Use 1 ml of the reagent per 1X 10⁶ of animal cells.

2.1.1.4 Blood: Mix 0.2 mL of blood and 0.8 mL RNA Extraction, vortex thoroughly.

2.1.2 PHASE SEPERATION

Centrifuge the samples at 12000rpm for 10 min at 4°C to remove cell debris. Transfer the supernatant to a new tube. Add 0.25 ml of chloroform per 1 ml of RNA Extraction. Cap sample tubes securely. Vortex samples vigorously for 15 seconds and incubate them at room temperature for 2 to 3 minutes. Centrifuge the samples at 12,000rpm for 10 minutes at 4°C. Following centrifugation, the mixture separates into lower red, phenol-chloroform phase, an interphase, and a colorless upper aqueous phase. RNA remains exclusively in the aqueous phase. Transfer upper aqueous phase carefully without disturbing the interphase into fresh tube.

2.1.3 RNA PRECIPITATION

Precipitate the RNA from the aqueous phase by mixing with isopropyl alcohol. Use 0.5 ml of isopropyl alcohol per 1 ml of RNA Extraction used for the initial homogenization. Incubate samples at -20°C for 15 minutes and centrifuge at 12,000rpm for 10 minutes at 4°C. The RNA precipitate, often invisible before centrifugation, forms a gel-like pellet on the side and bottom of the tube.

2.1.4 RNA WASH

Remove the supernatant completely. Wash the RNA pellet once with 75% ethanol, adding at least 1 ml of 75% ethanol per 1 ml of RNA Extraction used for the initial homogenization. Mix the samples by vortexing and centrifuge at 12,000rpm for 10 minutes at 4 °C. Remove all leftover ethanol.

2.1.5 REDISSOLVING RNA

Air-dry RNA pellet for 5-10 minutes. It is important not to let the RNA pellet dry completely as this will greatly decrease its solubility. Dissolve RNA in DEPC-treated water by passing solution a few times through a pipette tip.

2.1.6 SPECTROPHOTOMETRIC ANALYSIS

Dilute 1 μ l of RNA with 39 μ l of DEPC-treated water (1:40 dilution). Using 2.5 μ l RNA, take OD at 260 nm and 280 nm to determine sample concentration and purity. The A260/A280 ratio should be above 1.8. Apply the convention that 1 OD at 260 equals 40 μ g/ml RNA.

2.2 First Strand cDNA Synthesis

After thawing, mix and briefly centrifuge the components of the kit. Store on ice.

2.2.1 Add the following reagents into a sterile, nuclease-free tube on ice in the indicated order: Template RNA 2 μ g, Primer Oligo (dT) 18 primer 0.5 μ L and Random Hexamer primer 0.5 μ L or gene-specific primer 1 μ L, 5X Reaction Buffer 4 μ L, Wanwu®RT Enzyme Mix^a 1 μ L, Water, nuclease-free to 20 μ L, Total volume 20 μ L

2.2.2 Optional. If the RNA template is GC-rich or contains secondary structures, mix gently, centrifuge briefly and incubate at 65°C for 5 min. Chill on ice, spin down and place the vial back on ice.

2.2.3 Mix gently and centrifuge briefly.

2.2.4 Incubate for 5 min at 25°C, 30 min at 42°C.

Note. For GC-rich RNA templates the reaction temperature can be increased up to 50°C.

2.2.5 Terminate the reaction by heating at 85°C for 5 seconds.

2.3 Preparation of PCR Master Mix

2.3.1 For each 15 μ L reaction, prepare the following reaction mix:

2×SYBR Green qPCR Master Mix (High ROX)	7.5 μ L
Primer(2.5 μ M)	1.5 μ L
cDNA template	2.0 μ L
ddH ₂ O	4.0 μ L

2.3.2 PCR amplification

Pre-denaturation	95°C, 10min
cycle (40 次)	95°C, 15s→60°C, 60s
Melt Curve	60°C→95°C, 0.3 °C/15s

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Hefei WANWU technology CO.,LTD

A = CT (target gene, sample) - CT (internal standard gene, sample)

B = CT (target gene, control) - CT (internal standard gene, control)

K = A-B

RNA Expression = 2^{-K}