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## Genomic DNA extraction Lab Report

### 1 Laboratory Equipment and Reagents

#### 1.1 Laboratory equipment

Equipment	Manufacturers	Model
Centrifuge	DragonLab	D3024R
Clean bench	Suzhou Antai Air Tech Co., Ltd	SW-CJ-1FD
Standard reagent type pure water meter	Qingdao Fulum Technology Co., Ltd	FBZ2001-up-p
Ultramicro spectrophotometer	Thermo	NanoDrop2000

#### 1.2 Reagents

Reagents	Manufacturers	Order
TIANamp Genomic DNA Kit	Tiagen Biotech(Beijing) CO.,LTD	DP304
Isopropyl alcohol	Sinopharm Group Chemical Reagent Co., Ltd.	80109218
Centrifuge tube	Wanwu	
TIP	Wanwu	

### 2 Genomic DNA Extraction Experimental steps

#### 2.1 Samples preparation

2.1.1 For blood, please use 200  $\mu$ l fresh, frozen or anticoagulant adding blood. If less than 200  $\mu$ l, please make up with buffer GA to 200  $\mu$ l.

2.1.2 If the sample is blood from poultry, birds, amphibians, of which red blood cells have nucleolus, the amount should be reduced to 5-20  $\mu$ l and adjust the volume to 200  $\mu$ l with buffer GA.

2.1.3 The adherent cells should be treated to cell suspension first, then centrifuge the cells for 1 min at 10,000 rpm ( $\sim 11,200 \times g$ ), then discard the flow-through and re-suspend cell pellet in 200  $\mu$ l buffer GA.

2.1.4 Animal tissue (spleen  $< 10$ mg) should be treated to cell suspension first, then centrifuge the cells for 1 min at 10,000 rpm ( $\sim 11,200 \times g$ ), then discard the flow-through and re-suspend cell pellet in 200  $\mu$ l buffer GA.

2.2 Add 20  $\mu$ l Proteinase K, mix thoroughly by vortex. If the sample is tissue: incubate at 56°C until the tissue is completely lysed.

2.3 Add 200  $\mu$ l Buffer GB to the sample, mix thoroughly by vortex, and incubate at 70°C for 10 min to yield a homogeneous solution. Briefly centrifuge the 1.5 ml microcentrifuge tube to remove drops from the inside of the lid.

2.4 Add 200  $\mu$ l ethanol (96-100%) to the sample, and mix thoroughly by vortex for 15 s. A white

precipitate may form on addition of ethanol. Briefly centrifuge the 1.5 ml microcentrifuge tube to remove drops from the inside of the lid.

2.5 Pipet the mixture from step 4 into the Spin Column CB3 (in a 2 ml collection tube) and centrifuge at 12,000 rpm ( $\sim 13,400 \times g$ ) for 30 s. Discard flow-through and place the spin column into the collection tube.

2.6 Add 500  $\mu$ l Buffer GD (Ensure ethanol (96-100%) has been added) to Spin Column CB3, and centrifuge at 12,000 rpm ( $\sim 13,400 \times g$ ) for 30 s, then discard the flow-through and place the spin column into the collection tube.

2.7 Add 600  $\mu$ l Buffer PW (Ensure ethanol (96-100%) has been added) to Spin Column CB3, and centrifuge at 12,000 rpm ( $\sim 13,400 \times g$ ) for 30 s. Discard the flow-through and place the spin column into the collection tube.

2.8 Repeat Step 7.

2.9 Centrifuge at 12,000 rpm ( $\sim 13,400 \times g$ ) for 2 min to dry the membrane completely. Note: The residual ethanol of buffer PW may have some affection in downstream application. 10. Place the Spin Column CB3 in a new clean 1.5 ml microcentrifuge tube, and pipet 50-200  $\mu$ l Buffer TE directly to the center of the membrane. Incubate at room temperature (15-25°C) for 2-5 min, and then centrifuge for 2 min at 12,000 rpm ( $\sim 13,400 \times g$ ).