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Wuhan servicebio technology CO.,LTD

## PCR amplification and Agarose Gel Electrophoresis

### 1 Laboratory equipment and reagents

#### 1.1 Laboratory equipment

Equipment	Manufacturers	Model
Centrifuge	DragonLab	D3024R
PCR instrument	Beijing Dongsheng Innovation Biotechnology Co., Ltd.	ETC811
Constant temperature oscillation shaker	Changzhou Aohua Instrument Co., Ltd	SHZ-82A
Clean bench	Suzhou Antai Air Tech Co., Ltd	SW-CJ-1FD
Electrophoresis	Servicebio	FW-600
Gel imaging system	Tanon Science & Technology Co., Ltd.	Tanon-1600R
Standard reagent type pure water meter	Qingdao Fulum Technology Co., Ltd	FBZ2001-up-p
Ultramicro spectrophotometer	Thermo	NanoDrop2000

#### 1.2 Reagents

Reagents	Manufacturers	Order
TIANamp Genomic DNA Kit	Tiangen Biotech(Beijing) CO.,LTD	DP304
2×Fast Pfu PCR Master Mix	Servicebio	G3305
Agarose	Servicebio	G5056
50*TAE	Servicebio	G3001
MarkerI DNA ladder	Servicebio	WGM1100
Anhydrous ethanol	Sinopharm Group Chemical Reagent Co., Ltd.	10009218
Centrifuge tube	Servicebio	
TIP	Servicebio	
Primer	Servicebio	

### 2 Experimental steps

#### 2.1 Genomic DNA Extraction Experimental steps

##### 2.1.1 Samples preparation

2.1.1.1 For blood, please use 200  $\mu$ l fresh, frozen or anticoagulant adding blood. If less than 200  $\mu$ l, please make up with buffer GA to 200  $\mu$ l.

2.1.1.2 If the sample is blood from poultry, birds, amphibians, of which red blood cells have nucleolus, the amount should be reduced to 5-20  $\mu$ l and adjust the volume to 200  $\mu$ l with buffer

GA.

2.1.1.3 The adherent cells should be treated to cell suspension first, then centrifuge the cells for 1 min at 10,000 rpm ( $\sim 11,200 \times g$ ), then discard the flow-through and re-suspend cell pellet in 200  $\mu$ l buffer GA.

2.1.1.4 Animal tissue (spleen  $< 10\text{mg}$ ) should be treated to cell suspension first, then centrifuge the cells for 1 min at 10,000 rpm ( $\sim 11,200 \times g$ ), then discard the flow-through and re-suspend cell pellet in 200  $\mu$ l buffer GA.

2.1.2 Add 20  $\mu$ l Proteinase K, mix thoroughly by vortex. If the sample is tissue: incubate at  $56^\circ\text{C}$  until the tissue is completely lysed.

2.1.3 Add 200  $\mu$ l Buffer GB to the sample, mix thoroughly by vortex, and incubate at  $70^\circ\text{C}$  for 10 min to yield a homogeneous solution. Briefly centrifuge the 1.5 ml microcentrifuge tube to remove drops from the inside of the lid.

2.1.4 Add 200  $\mu$ l ethanol (96-100%) to the sample, and mix thoroughly by vortex for 15 s. A white precipitate may form on addition of ethanol. Briefly centrifuge the 1.5 ml microcentrifuge tube to remove drops from the inside of the lid.

2.1.5 Pipet the mixture from step 4 into the Spin Column CB3 (in a 2 ml collection tube) and centrifuge at 12,000 rpm ( $\sim 13,400 \times g$ ) for 30 s. Discard flow-through and place the spin column into the collection tube.

2.1.6 Add 500  $\mu$ l Buffer GD (Ensure ethanol (96-100%) has been added) to Spin Column CB3, and centrifuge at 12,000 rpm ( $\sim 13,400 \times g$ ) for 30 s, then discard the flow-through and place the spin column into the collection tube.

2.1.7 Add 600  $\mu$ l Buffer PW (Ensure ethanol (96-100%) has been added) to Spin Column CB3, and centrifuge at 12,000 rpm ( $\sim 13,400 \times g$ ) for 30 s. Discard the flow-through and place the spin column into the collection tube.

2.1.8 Repeat Step 7.

2.1.9 Centrifuge at 12,000 rpm ( $\sim 13,400 \times g$ ) for 2 min to dry the membrane completely. Note: The residual ethanol of buffer PW may have some affection in downstream application

2.1.10 Place the Spin Column CB3 in a new clean 1.5 ml microcentrifuge tube, and pipet 50-200  $\mu$ l Buffer TE directly to the center of the membrane. Incubate at room temperature ( $15\text{-}25^\circ\text{C}$ ) for 2-5 min, and then centrifuge for 2 min at 12,000 rpm ( $\sim 13,400 \times g$ ).

## 2.2 PCR amplification

2.2.1 0.2ml PCR tube, 20ul pcr amplification system.

2 x Fast Pfu PCR Master Mix	25 $\mu$ l
Forward Primer (10 $\mu$ M)	1 $\mu$ l
Reverse Primer (10 $\mu$ M)	1 $\mu$ l
Template	Variable
ddH <sub>2</sub> O	Add to 50 $\mu$ L

## 2.2.2 PCR amplification

Pre-denaturation	98°C, 2min	} 30×cycle
Denatured	98°C, 20s	
Annealing	55°C, 20s	
Extended	72°C, 10s	
End extension	72°C, 5min	

## 3 Agarose gel electrophoresis:

3.1 Melt 3g agarose in 100ml of 1X TAE in the microwave until all agarose is dissolved and there are no stringy pieces. Cool the liquid under cold running water for 10-15 seconds or allow to cool by letting it sit at RT until it is not too hot to hold. Add a very small amount of ethidium bromide (1ul per 20ml is sufficient). You can do this simply by dipping the pipet tip into the stock solution and then swirling this into the liquid agar. Ethidium bromide is used to visualize the DNA when viewed under UV light.

3.2 Pour the gel into the gel mold held in place by the clamp, with the desired comb in place. Immediately rinse the flask in RO water and place on drying rack. Allow the gel to dry for about 15 to 20 minutes.

3.3 Remove comb and transfer gel in gel mold to gel box with TBE buffer, making sure that the gel is completely submerged (do not fill past max fill line).

3.4 Load the gel. Each sample needs loading buffer, which should be used at a ratio of 1:6. Loading buffer is used as a marker for migration rates and also to give the sample density so that it sinks to the bottom of the well.

3.5 After all samples and ladder(s) have been loaded, place the lid on the gel box, make sure it is plugged into the powersource, and turn the powersource on. Run it at about 92 volts for 30-45 minutes, or until loading buffer bands are where you want them.

3.6 Check the gel when it is done using the Gel imaging system. If it's not spread out enough, run longer. Remember to turn off the power when you disconnect the lid from the box! Also, when not in use, please make sure the lid is on to minimize evaporation of buffer and to keep dust out. When you have taken a picture or saved the image, throw the gel away in the trash.