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Immunohistochemical Experiment Report of Paraffin Section

I Experimental equipments and reagents

1. Experimental equipments

| Name | Manufacturer | Model |
|-----------------------|---|---------------------|
| Dehydrator | DIAPATH | Donatello |
| Embedding machine | Wuhan Junjie Electronics Co., LTD | JB-P5 |
| Pathological slicer | Shanghai Leica Instrument Co., LTD | RM2016 |
| Refrigerating table | Wuhan Junjie Electronics Co., LTD | JB-L5 |
| Organizing spreader | Zhejiang Jinhua Kedi Instrument Equipment Co., LTD | KD-P |
| Oven | Shanghai Huitai Instrument Manufacturing Co., LTD | DHG-9140A |
| Slide | Wanwu | |
| Cover glass | Jiangsu Shitai Experimental Equipment Co., LTD | 10212432C |
| Microwave oven | Galanz Microwave Oven Electrical | P70D20TL-P4 |
| Decolorization shaker | Appliances co., LTD | TSY-B |
| Vortex mixer | Wanwu | MX-F |
| Palm centrifuge | Wanwu | D1008E |
| Pipette gun | Dragon | KE0003087/KA0056573 |
| Tissue pencil | Wanwu | WG1066-1 |
| Microscope | Nikon | E100 |
| Imaging system | Japan Nikon | Nikon DS-U3 |

2. Main experimental reagents

| Reagent | Manufacturer | Model | Dilution Ratio |
|-------------------|---|-----------|----------------|
| Anhydrous ethanol | China National Pharmaceutical Group Chemical Reagent Co., LTD | 100092683 | |

| | | |
|---|---|-----------|
| Xylene | Shanghai Lingfeng Chemical Reagent Co., LTD | 1330-20-7 |
| Citric acid(PH6.0)antigen repair solution | Wanwu | G1202 |
| Phosphate buffer saline | Wanwu | G0002 |
| 4% Paraformaldehyde | Wanwu | G1101 |
| 3% Hydrogen peroxide | China National Pharmaceutical Group Chemical Reagent Co., LTD | 10011208 |
| Bovine serum albumin | Wanwu | G5001 |
| Normal rabbit serum | Wanwu | G1209 |
| Hematoxylin stain solution | Wanwu | G1004 |
| Hematoxylin differentiation solution | Wanwu | G1309 |
| Hematoxylin returning blue solution | Wanwu | G1340 |
| Glycerine gelatin sealant | Wanwu | G1402 |
| Primary antibody: | | |
| Secondary antibody: | | |
| Immunohistochemical kit | Wanwu | G1211 |
| AEC chromogenic agent | | |

II Immunohistochemical experimental procedure of paraffin section

1. Deparaffinzing and rehydrating the paraffin section: put the sections into xylene I for 15 minutes--xylene II for 15 minutes--xylene III for 15 minutes--absolute ethanol I for 5 minutes--absolute ethanol II for 5 minutes--85% alcohol for 5 minutes--75% alcohol for 5 minutes--rinse in distilled water.
2. Antigen retrieval: The tissue sections are placed in a repair box filled with citric acid(PH6.0) antigen retrieval buffer for antigen retrieval in a microwave oven, heated on medium power for 8 minutes until boiling, then turned off the microwave oven, kept warm for 8 minutes and then transferred to medium-low power for heating 7minutes. During this process, excessive evaporation of buffer should be prevented and the sections should not be allowed to dry. To cool to room temperature before proceeding, the sections are placed in PBS(PH7.4) and shaken on the decolorization shaker 3 times for 5 minutes each.
3. Blocking endogenous peroxidase activity: the sections are placed in 3% hydrogen peroxide and incubated at room temperature in darkness for 25 minutes. The sections are placed in

-
- PBS(PH7.4) and shaken on a decolorizing shaper 3 times for 5 minutes each.
4. Serum sealing: 3%BSA was added to the circle to evenly cover the tissue, and the tissues are sealed for 30 minutes at room temperature. (Primary antibody is sealed with normal rabbit serum from goat source and other sources are sealed with BSA).
 5. Primary antibody incubation: the sealing solution is gently removed, the primary antibody prepared with PBS(PH7.4) in a certain proportion is added to the sections, and the sections are placed flat in a wet box and incubated overnight at 4°C. (Add a small amount of water in the wet box to prevent evaporation of antibodies).
 6. Secondary antibody incubation: the sections are placed in PBS(PH7.4) and washed by shaking on the decolorizing shaker 3 times for 5 minutes each. After the sections are slightly shaken and dried, the tissues are covered with secondary antibody (HRP labeled) from the corresponding species of primary antibody and incubated at room temperature for 50 minutes.
 7. AEC chromogenic reaction: the sections are placed in PBS(PH7.4) and shaken on the decoloring shaker 3 times for 5 minutes each. ACE color developing solution newly prepared is added in the circle after the sections are slightly dried. The color developing time is controlled under the microscope. The positive is brownish yellow. Rinse the sections with tap water to stop the reaction.
 8. Nucleus counterstaining: the sections are counterstained with hematoxylin stain solution for about 3 minutes; washed with tap water; differentiated with hematoxylin differentiation solution for several seconds; washed with tap water; treated with hematoxylin returning blue solution; washed with running water.
 9. Stabilizing with mounting medium: after the nuclear microscopy is completed, the section is placed in pure water, slightly dried, then mount the sections with neutral gum and air-dried naturally.
 10. Visualize staining of tissue under a microscope, acquisitive and analysis image.

III Interpretation of the immunohistochemical results of paraffin section

The nucleus of hematoxylin stained is blue, and the positive expression of AEC is red.