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Immunohistochemical Experiment Report of Cell Drop Tablet

I Experimental equipment and reagents

1.Experimental equipment

Name	Manufacturer	Model
Slide	Wanwu	Slide
Cover glass	Jiangsu Shitai Experimental Equipment Co., LTD.	10212432C
Vortex mixer	Wanwu	MX-F
Palm centrifuge	Wanwu	D1008E
Decolorization shaker	Wanwu	TSY-B
Pipette gun	Dragon	KE0003087/KA0056573
Tissue pencil	Wanwu	WG1066-1
Microscope	Nikon	E100
Imaging system	Japan Nikon	Nikon DS-U3

2.Main experimental reagents

Reagent	Manufacturer	Model	Dilution Ratio
Anhydrous ethanol	China National Pharmaceutical Group Chemical Reagent Co., LTD	100092683	
Xylene	Shanghai Lingfeng Chemical Reagent Co., LTD	1330-20-7	
Phosphate buffer saline	Wanwu	G0002	
Bovine serum albumin	Wanwu	G5001	
Rupture fluid	Wanwu	G1204	
4% Paraformaldehyde	China National Pharmaceutical	G1101	
Neutral gum	Group Chemical Reagent Co., LTD	10004160	
Primary antibody:			
Secondary antibody:			

Immunohistochemical kit

Wanwu

G1211

DAB chromogenic agent

II Experimental procedure of immunocytochemistry

1. The cell suspension is centrifuged at 2800rpm at 4°C for 5 minutes, the supernatant is discarded, and 2mL 4% Paraformaldehyde is added for fixation according to the amount of cells deposited at the bottom. If the cell pellet is not visible to the naked eye, it is centrifuged at 4°C for 10 minutes at 3000rpm.
2. Take the cell suspension resuspended in 4% Paraformaldehyde, centrifuge at 2800rpm at 25°C for 5 minutes, abandon the supernatant, according to the bottom of the precipitation to add PBS(PH7.4): ①If the cell pellet is not visible to the naked eye, it is centrifuged at 25°C for 10 minutes at 3000rpm., when there is still no precipitation, leave about 0.2mL of liquid at the bottom, add 0.5mL of PBS and mix well, then use a pipette gun to suck up 200μL, and drop it into a small circle (3cm×2cm ellipse) drawn in advance with a tissue pencil; ②if the amount of cell precipitation is very small and the similar in size to mung beans, discard the supernatant, keep the bottom precipitation, add 0.5mL PBS, blow it with a pipette gun and mix it well, then absorb 200μL, and spread it in the small circle (3cm×2cm ellipse) drawn in advance with the composition; ③if the amount of cell precipitation is large and the size of soybean is similar or larger, the supernatant should be discarded and the bottom precipitation should be retained, add 2mL PBS, mix with a pipette gun, draw 150μL, and apply it in advance with a group pen in a large circle(maximum about 5cm long side, short edge of about 2cm of the elliptic), how much the amount of cells under a microscope, if it is still too thick, dilute the cell suspension double-fold with PBS until the amount of cells is evenly spread under the microscope
3. The cell suspension with pipette gun covered the whole circle, and place the section to dry naturally.
4. Cell fixation and membrane rupture: draw a circle on the position where the cells are smeared on the slide with a chemical pen (to prevent the antibody from flowing away), and add 50-100μL of 4% Paraformaldehyde. After fixation at room temperature for 20 minutes, wash away the fixative, add Rupture working fluid, incubate at room temperature for 20 minutes, and wash 3 times with PBS for 5 minutes each time.
5. Serum sealing: 3%BSA was added to the circle to evenly cover the tissue, and the tissues are sealed for 30 minutes at room temperature. (Primary antibody is sealed with normal rabbit serum from goat source and other sources are sealed with BSA).
6. Primary antibody incubation: the sealing solution is gently removed, the primary antibody prepared with PBS(PH7.4) in a certain proportion is added to the sections, and the sections are placed flat in a wet box and incubated overnight at 4°C.

7. Secondary antibody incubation: the sections are placed in PBS(PH7.4) and washed by shaking on the decolorizing shaker 3 times for 5 minutes each. After the sections are slightly shaken and dried, the tissues are covered with secondary antibody from the corresponding species of primary antibody and incubated at room temperature for 50 minutes.
8. DAB chromogenic reaction: the sections are placed in PBS(PH7.4) and shaken on the decoloring shaker 3 times for 5 minutes each. DAB color developing solution newly prepared is added in the circle after the sections are slightly dried. The color developing time is controlled under the microscope. The positive is brownish yellow. Rinse the sections with tap water to stop the reaction.
9. Nucleus counterstaining: the sections are counterstained with hematoxylin stain solution for about 3 minutes; washed with tap water; differentiated with hematoxylin differentiation solution for several seconds; washed with tap water; treated with hematoxylin returning blue solution; washed with running water.
10. Dehydration and mounting: place the section in 75% alcohol for 5 minutes--85% alcohol for 5 minutes--absolute ethanol I for 5 minutes--anhydrous ethanol II for 5 minutes--xylene I for 5 minutes, dehydrated and transparent, remove the sections from xylene and let them dry slightly, then mount the sections with neutral gum.
11. Visualize staining of tissue under a microscope, acquisitive and analysis image.

III Interpretation of the immunocytochemistry results

The nucleus of hematoxylin stained is blue, and the positive expression of DAB is brownish yellow.