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1.2 Major reagents

Name	Producer	Code	Dilution
Ethanol	SCRC	100092183	

Xylene	上海凌峰化学试剂有限公司	1330-20-7
EDTA antigen retrieval solution (pH 8.0)	Wanwu	G1206
EDTA antigen retrieval solution (pH 9.0)	Wanwu	G1203
Sodium citrate antigen retrieval solution (pH 6.0)	Wanwu	G1202
PBS solution	Wanwu	G0002
3%H ₂ O ₂	SCRC	
FITC-TSA	Wanwu	G1222
CY3-TSA	Wanwu	G1223
647-TSA	Wanwu	G1224
Spontaneous fluorescence quenching reagent	Wanwu	G1221
BSA primary antibody:	Wanwu	G5001
secondary antibody		
DAPI	Wanwu	G1012

2 Procedure

2.1 Deparaffinize and rehydrate: incubate sections in 2 changes of xylene, 15 min each. Dehydrate in 2 changes of pure ethanol for 5 min, followed by dehydrate in gradient ethanol of 85% and 75% ethanol, respectively, 5 min each. Wash in distilled water.

2.2 Antigen retrieval: immerse the slides in EDTA antigen retrieval buffer (pH 8.0) and maintain at a sub-boiling temperature for 8 min, standing for 8 min and then followed by another sub-boiling temperature for 7 min. Be sure to prevent buffer solution evaporate. Let air cooling. Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Use the right antigen retrieval buffer and heat extent according to tissue characteristics.

2.3 Circle, Block endogenous peroxidase : wash three times with PBS (pH 7.4) in a Rocker device, 5 min each eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Immerse in 3% H₂O₂ and incubate at room temperature for 15 min, kept in dark place. Then wash again three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.4 Block with serum: eliminate obvious liquid, mark the objective tissue with liquid blocker pen.

Cover objective tissues with 10% donkey serum (for the case of primary antibody originated from goat) or 3% BSA (for the case of primary antibody originated from others) at room temperature for 30 min.

2.5 First primary antibody :throw away the blocking solution slightly. Incubate slides with the first primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

2.6 Corresponding secondary antibody marked with HRP:wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to first primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.7 CY3-TSA solution:wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each.Incubate slides with CY3-TSA solution(diluted with TBST appropriately) for 10 min in dark condition.After that,wash slides three times with TBST in a Rocker device, 5 min each.

2.8 Microwave treatment:immerse the slides in EDTA antigen retrieval buffer (pH 8.0) and maintain at a sub-boiling temperature for 8 min, standing for 8 min and then followed by another sub-boiling temperature for 7 min,to remove the primary antibodies and Secondary antibodies combined with tissue.Be sure to prevent buffer solution evaporate.

2.9 Second primary antibody :Incubate slides with primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

2.10 Second corresponding secondary antibody marked with HRP :wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to second primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.11 FITC-TSA solution:wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each.Incubate slides with FITC-TSA solution(diluted with TBST appropriately) for 10 min in dark condition.After that,wash slides three times with TBST in a Rocker device, 5 min each.

2.12 Microwave treatment:immerse the slides in EDTA antigen retrieval buffer (pH 8.0) and maintain at a sub-boiling temperature for 8 min, standing for 8 min and then followed by another sub-boiling temperature for 7 min,to remove the primary antibodies and Secondary antibodies combined with tissue.Be sure to prevent buffer solution evaporate.

2.13 Incubate the third primary antibody :Incubate slides with primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

2.14 The third corresponding secondary antibody marked with HRP :wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to second primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.15 647-TSA solution:wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min

each. Incubate slides with 647-TSA solution (diluted with TBST appropriately) for 10 min in dark condition. After that, wash slides three times with TBST in a Rocker device, 5 min each.

2.16 Microwave treatment: immerse the slides in EDTA antigen retrieval buffer (pH 8.0) and maintain at a sub-boiling temperature for 8 min, standing for 8 min and then followed by another sub-boiling temperature for 7 min, to remove the primary antibodies and Secondary antibodies combined with tissue. Be sure to prevent buffer solution evaporate.

2.17 Incubate the fourth primary antibody : Incubate slides with primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

2.18 The fourth corresponding secondary antibody marked with 594 labeled Fluorescence secondary antibody : wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to second primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.19 DAPI counterstain in nucleus: incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.20 Spontaneous fluorescence quenching: eliminate obvious liquid, incubate slides with spontaneous fluorescence quenching reagent for 5 min, then wash slides under flowing water for 10 min.

2.21 Mount: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Throw away liquid slightly, then coverslip with anti-fade mounting medium.

2.22 Detection and collect images by slice scanner : DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm. 647 glows pink by excitation wavelength 608-648 nm and emission wavelength 672-712 nm. (647 was originally red, in order to distinguish it from CY3, we set it to pink light.). 594 glows purplish red by excitation wavelength 594 nm and emission wavelength 615 nm. (594 was originally red, in order to distinguish it from CY3, we set it to purplish red light.)

3 Results

Nucleus is blue by labeling with DAPI. Positive cells are green, red, pink or purplish red according to the fluorescent labels used.