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Immunofluorescence protocol (paraffin-slides)**1. Apparatus and reagents**

1.1 Apparatus

Name	Producer	Model
Dehydrator	DIAPATH	Donatello
Paraffin embedding machine	武汉俊杰电子有限公司	JB-P5
pathology slicer	Leica	RM2016
Frozen machine	武汉俊杰电子有限公司	JB-L5
Water bath-slide drier	浙江省金华市科迪仪器设备有限公司	KD-P
Oven	上海慧泰仪器制造有限公司	DHG-9140A
Glass microscope slides	Wanwu	
Coverslips	CITOTEST	10212432C
Microwave	Glanze	P70D20TL-P4
Rocker	Wanwu	TSY-B
Vortex	Wanwu	MX-F
Micro-centrifuge	Wanwu	D1008E
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Wanwu	WG1066-1
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

1.2 Major reagents

Name	Producer	Code	Dilution
Ethanol	SCRC	100092183	

Xylene	上海凌峰化学试剂有限公司	1330-20-7
EDTA antigen retrieval solution (pH 8.0)	Wanwu	G1206
EDTA antigen retrieval solution (pH 9.0)	Wanwu	G1203
Sodium citrate antigen retrieval solution (pH 6.0)	Wanwu	G1202
PBS solution	Wanwu	G0002
3%H ₂ O ₂	SCRC	
TSA-FITC	Wanwu	G1222
TSA-CY3	Wanwu	G1223
Spontaneous fluorescence quenching reagent	Wanwu	G1221
BSA primary antibody:	Wanwu	G5001
secondary antibody		
DAPI	Wanwu	G1012
anti-fade mounting medium	Wanwu	G1401

2 Procedure

2.1 Deparaffinize and rehydrate: incubate sections in 2 changes of xylene, 15 min each. Dehydrate in 2 changes of pure ethanol for 5 min, followed by dehydrate in gradient ethanol of 85% and 75% ethanol, respectively, 5 min each. Wash in distilled water.

2.2 Antigen retrieval: immerse the slides in EDTA antigen retrieval buffer (pH 8.0) and maintain at a sub-boiling temperature for 8 min, standing for 8 min and then followed by another sub-boiling temperature for 7 min. Be sure to prevent buffer solution evaporate. Let air cooling. Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Use the right antigen retrieval buffer and heat extent according to tissue characteristics.

2.3 Circle, Block endogenous peroxidase : wash three times with PBS (pH 7.4) in a Rocker device, 5 min each eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Immerse in 3% H₂O₂ and incubate at room temperature for 25 min, kept in dark place. Then wash again three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.4 Block with serum: eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Cover objective tissues with 10% donkey serum (for the case of primary antibody originated from goat) or 3% BSA (for the case of primary antibody originated from others) at room temperature for 30 min.

2.5 First primary antibody :throw away the blocking solution slightly. Incubate slides with the first primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

2.6 Corresponding secondary antibody marked with HRP:wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to first primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.7 TSA-FITC(or TSA-CY3) solution:wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each.Incubate slides with TSA-FITC solution(diluted with TBST appropriately) for 10 min in dark condition.After that,wash slides three times with TBST in a Rocker device, 5 min each.

2.8 Microwave treatment:immerse the slides in EDTA antigen retrieval buffer (pH 8.0) and maintain at a sub-boiling temperature for 8 min, standing for 8 min and then followed by another sub-boiling temperature for 7 min,to remove the primary antibodies and Secondary antibodies combined with tissue.Be sure to prevent buffer solution evaporate.

2.9 Second primary antibody :Incubate slides with primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

2.10 Second corresponding secondary antibody :wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to second primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.11 Spontaneous fluorescence quenching:eliminate obvious liquid,incubate slides with spontaneous fluorescence quenching reagent for 5 min,then wash slides under flowing water for 10 min.

2.12 DAPI counterstain in nucleus:incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.13 Mount: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Throw away liquid slightly, then coverslip with anti-fade mounting medium.

2.14 Microscopy detection and collect images by Fluorescent Microscopy. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

3 Results

Nucleus is blue by labeling with DAPI. Positive cells are green or red according to the fluorescent labels used.