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## Phalloidine staining protocol for cell sample (FITC)

### 1. Apparatus and reagents

#### 1.1 Apparatus

Name	Producer	Model
Glass microscope slides	Wanwu	
Vortex	Wanwu	MX-F
Micro-centrifuge	Wanwu	D1008E
Rocker	Wanwu	TSY-B
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Wanwu	WG1066-1
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

#### 1.2 Major reagents

Name	Producer	Code	Dilution
PBS solution	Wanwu	G0002	
Permeabilize solution	Wanwu	G1204	
phalloidin (FITC)	Wanwu	G1028	1:200
DAPI	Wanwu	G1012	
anti-fade mounting medium	Wanwu	G1401	

### 2 Procedure

2.1 Permeabilization: dry the cell climbing slides slightly. Mark the objective area with liquid blocker pen, where add 50-100  $\mu$ l of permeabilize working solution. Incubate for 20 min at room temperature. Wash three times with PBS solution, 5 min each.

2.2 Phalloidine (green) stain: add phalloidine diluted solution to cover with the objective area, incubate at room temperature for 2 h in the dark.

2.3 DAPI counterstain in nucleus: wash three times with PBS (pH 7.4) in a Rocker device, 5 min

each. Then incubate with DAPI solution at room temperature for 10 min in the dark.

2.4 Mount: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Throw away liquid slightly. Put the slides on a glass microscope slide and then mount with anti-fade mounting medium

2.5 Microscopy detection and collect images by Fluorescent Microscopy. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

### **3 Results**

Nucleus is blue by labeling with DAPI. Positive cells labelled by Phalloidine (FITC) are green