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Immunofluorescence protocol (cell climbing slides)

1. Apparatus and reagents

1.1 Apparatus

Name	Producer	Model
Glass microscope slides	Wanwu	
Vortex	Wanwu	MX-F
Micro-centrifuge	Wanwu	D1008E
Rocker	Wanwu	TSY-B
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Wanwu	WG1066-1
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

1.2 Major reagents

Name	Producer	Code	Dilution
PBS solution	Wanwu	G0002	
Permeabilize solution	Wanwu	G1204	
BSA	Wanwu	G5001	
Primary antibody			
Secondary antibody			
DAPI	Wanwu		
anti-fade mounting medium	Wanwu	G1401	

2 Procedure

2.1 Permeabilization: dry the cell climbing slides slightly. Mark the objective area with liquid blocker pen, where add 50-100 μ l of permeabilize working solution. Incubate for 20 min at room temperature. Wash three times with PBS solution, 5 min each.

2.2 Block with serum: eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Cover objective area with 10% donkey serum (for the case of primary antibody originated from goat) or 3% BSA (for the case of primary antibody originated from others) at room temperature

for 30 min.

2.3 Primary antibody: remove the blocking solution. Incubate cells with primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box.

2.4 Secondary antibody: wash cell climbing slides three times with PBS (pH 7.4), placed in a Rocker device and 5 min each. Cover cell climbing slides with secondary antibody (appropriately respond to primary antibody in species) , incubate at room temperature for 50 min.

2.5 DAPI counterstain in nucleus: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.6 Mount: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Throw away liquid slightly. Put the slides on a glass microscope slide and then mount with anti-fade mounting medium.

2.7 Microscopy detection and collect images by Fluorescent Microscopy.: DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

3 Results

Nucleus is blue by labeling with DAPI. Positive cells are green or red according to the fluorescent labels used.