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Paraffin Thioflavin staining protocol

1. Apparatus and reagents

1.1 Apparatus

Name	Producer	Model
Dehydrator	DIAPATH	Donatello
Paraffin embedding machine	武汉俊杰电子有限公司	JB-P5
pathology slicer	Leica	RM2016
Frozen machine	武汉俊杰电子有限公司	JB-L5
Water bath-slide drier	浙江省金华市科迪仪器设备有限公司	KD-P
Oven	上海慧泰仪器制造有限公司	DHG-9140A
Glass microscope slides	Wanwu	
Coverslips	CITOTEST	10212432C
Rocker	Wanwu	TSY-B
Vortex	Wanwu	MX-F
Micro-centrifuge	Wanwu	D1008E
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Wanwu	WG1066-1
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

1.2 Major reagents

Name	Producer	Code	Dilution
Ethanol	SCRC	100092183	
Xylene	上海凌峰化学试剂有限公司	1330-20-7	
PBS solution	Wanwu	G0002	
DAPI	Wanwu	G1012	
anti-fade mounting medium	Wanwu	G1401	

2 Procedure

2.1 Deparaffinize and rehydrate: incubate sections in 2 changes of xylene, 15 min each. Dehydrate in 2 changes of pure ethanol for 5 min, followed by dehydrate in gradient ethanol of 85% and 75% ethanol, respectively, 5 min each. Wash in distilled water

2.2 Thioflavin staining: Prepared with 50% alcohol at a concentration of 0.3%, filtered and incubated at room temperature for 8min.

2.3 DAPI counterstain in nucleus: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.4 Mount: Throw away liquid slightly, then coverslip with anti-fade mounting medium.

2.5 Microscopy detection and collect images by Fluorescent Microscopy. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

3 Results

Nucleus is blue by labeling with DAPI. Thioflavin positive cells labelled by fluorescein are green.