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cell climbing slides edu staining protocol

1. Apparatus and reagents

1.1 Apparatus

Name	Producer	Model
Glass microscope slides	Wanwu	
Vortex	Wanwu	MX-F
Micro-centrifuge	Wanwu	D1008E
Rocker	Wanwu	TSY-B
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Wanwu	WG1066-1
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

1.2 Major reagents

Name	Producer	Code	Dilution
PBS solution	Wanwu	G0002	
Permeabilize solution	Wanwu	G1204	
Glycine	Wanwu	G5010-500G	
EDU staining kit	Wanwu	C10310	
DAPI	Wanwu	C10338	
anti-fade mounting medium	Wanwu	G1401	

2 Procedure

2.1 Add 2 mg/mL of glycine to each well, incubate in rocker device for 5 minutes, and discard the glycine solution.

Note: The action of this step is to neutralize paraformaldehyde, which can be omitted when Cell Fixation is performed in other ways.

2.2 Add 100 μ L of PBS to each well, wash in rocker device for 5 minutes, and discard PBS;

2.3 Add 100 μ L of penetrant to each well (0.5% TritonX-100 of PBS), incubate in rocker device for 10 minutes; wash with PBS once for 5 minutes.

Note: When the experiment requires other antibody staining, or because some cell types are highly adsorptive to dyes, it may be necessary to enhance cell membrane permeability.

2.4 Add 100 μ L 1X Apollo[®] staining reaction solution to each well, and discard the staining reaction solution after incubating for 30 minutes at room temperature and protected from light;

2.5 Add 100 μ L penetrant (0.5% TritonX-100 in PBS) and wash in a Rocker device 2 to 3 times, 10 minutes each time, discard the penetrant;

2.6 (Reinforced)add 100 μ L of methanol to each hole 1 or 2 times for 5 minutes each time; wash once with PBS for 5 min each time.

2.7 DAPI counterstain in nucleus: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.8 Mount: Throw away liquid slightly, then coverslip with anti-fade mounting medium

2.9 Microscopy detection and collect images by Fluorescent Microscopy. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

3 Results

Nucleus is blue by labeling with DAPI. Positive cells are green or red according to the fluorescent labels used.