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**DABTunel assay report (Frozen-slides)****1. Apparatus and reagents**

## 1.1 Apparatus

<b>Name</b>	<b>Producer</b>	<b>Model</b>
Freezing microtome	Thermo	Cryotome E
Glass microscope slides	wanwu	
Coverslips	CITOTEST	10212432C
Rocker	wanwu	TSY-B
Vortex	wanwu	MX-F
Micro-centrifuge	wanwu	D1008E
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	wanwu	WG1066-1
Refrigerator	Haier	BCD-192TGN
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

## 1.2 Major reagents

<b>Name</b>	<b>Producer</b>	<b>Model</b>
4% <b>paraformaldehyde</b>	wanwu	G1101
PBS solution	wanwu	G0002
Proteinase K	wanwu	G1205
Permeabilize solution	wanwu	G1204
3% H2O2	wanwu	G0115
Hematoxylin staining solution	wanwu	G1004
Hematoxylin solution	differentiate wanwu	G1039
Hematoxylin bluing solution	wanwu	G1040
Resin mounting medium	SCRC	<b>10004160</b>

---

Tunel Assay Kit	wanwu	G1507
DAB reagent	wanwu	G1212

---

## 2. Procedure

2.1 Fix frozen-slides: restore frozen-slides to room temperature and dry in air. Fix in cold acetone 10 min and then dry in air. Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.2 Antigen retrieval: eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Add proteinase K working solution to cover objectives and incubate at 37°C for 20 min. then wash three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.3 ( optional step ) Permeabilization: eliminate excess liquid, add permeabilize working solution to cover objective tissue, then incubate at room temperature for 20 min. wash three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.4 Block endogenous peroxidase: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Immerse in 3% H<sub>2</sub>O<sub>2</sub> and incubate at room temperature for 20 min, kept in dark place. Then wash again three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.5 Equilibrium at room temperature: After the slices are slightly dried, buffer is added to the tissues in the circle, and the buffer is incubated at room temperature for 10 minutes

2.6 Tunel reaction: Take appropriate amount of TDT enzyme, dUTP and buffer in the tunel kit according to the number of slices and tissue size and mix at 1:5:50 ratio。 Prepare this reaction solution according to demand before use. Add this mixture to objective tissue placed in a flat wet box, incubate at 37°C for 1 h. be sure to keep the wet box moist by adding water.

2.7 Reagent Streptavidin-HRP and TBST are mixed at a ratio of 1:200, added to the circle to cover the tissue, the slices are placed in a wet box, and incubated in a 37°C incubator for 30 minutes. The slides were placed in PBS (PH7.4) and washed with shaking on a decolorizing shaker 3 times, 5 min each time.

2.8 DAB developing: dry sections slightly, add fresh prepared DAB chromogenic reagent to marked tissue. Manage reaction time by observing in microscopy until nucleus shows brown-yellow. Then stop developing reaction by wash in pure water.

2.9 Counterstain in nucleus with Hematoxylin staining solution for 1 min and wash in pure water.

---

---

Treat with the differentiate solution for a few seconds, wash in pure water. Back to blue by bluing solution, wash in pure water.

2.10 Dehydrate the slices in 4 cylinders of 100% alcohol for 5 minutes each, then soak them in n-butanol for 5 minutes, Put it in xylene for transparency, more than 5min. Dry briefly and mount with resin mounting medium.

### 3. Results

Nucleus stained with hematoxylin are blue. The positive apoptosis cells developed by DAB reagent have brown-yellow nucleus.