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(FITC) tunel assay protocol (frozen-slides)

1. Apparatus and reagents

1.1 Apparatus

Name	Producer	Model
Freezing microtome	Thermo	Cryotome E
Glass microscope slides	wanwu	
Coverslips	CITOTEST	10212432C
Rocker	wanwu	TSY-B
Vortex	wanwu	MX-F
Micro-centrifuge	wanwu	D1008E
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	wanwu	WG1066-1
Refrigerator	Haier	BCD-192TGN
Ortho-Fluorescent Microscopy	Nikon	NIKON ECLIPSE C1
Imaging system	Nikon	NIKON DS-U3

1.2 Major reagents

Name	Producer	Model
4% paraformaldehyde	wanwu	G1101
PBS solution	wanwu	G0002
Proteinase K	wanwu	G1205
Permeabilize solution	wanwu	G1204
DAPI	wanwu	G1012
anti-fade mounting medium	wanwu	G1401
Tunel assay kit	wanwu	G1501

2. Procedure

2.1 Fix frozen-slides: restore frozen-slides to room temperature and dry in air. Fix in cold acetone 10 min and then dry in air. Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each

2.2 Antigen retrieval: eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Add proteinase K working solution to cover objectives and incubate at 37°C for 25 min. then wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. (Method for configuring working solution of proteinase K, stock solution: PBS=1:9)

2.3 Permeabilization: eliminate excess liquid, add permeabilize working solution to cover objective tissue, then incubate at room temperature for 20 min. wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. (The membrane breaking fluid is 0.1% triton. Configuration method, triton stock solution: PBS=1:1000)

2.4 Equilibrium at room temperature: After the slices are slightly dried, buffer is added to the tissues in the circle, and the buffer is incubated at room temperature for 10 minutes

2.5 TUNEL reaction: Take appropriate amount of TDT enzyme, dUTP and buffer in the TUNEL kit according to the number of slices and tissue size and mix at 1:5:50 ratio. Prepare this reaction solution according to demand before use. Add this mixture to objective tissue placed in a flat wet box, incubate at 37°C for 2 h. be sure to keep the wet box moist by adding water.

2.6 DAPI counterstain in nucleus: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min, kept in dark place

2.7 Mount: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Throw away liquid slightly, then coverslip with anti-fade mounting medium.

2.8 Microscopic examination and collecting images through fluorescence microscope. DAPI emits blue light at an ultraviolet excitation wavelength of 330-380 nm and an emission wavelength of 420 nm; FITC has an excitation wavelength of 465-495 nm and an emission wavelength of 515-555 nm, and emits green.

3. Results

Nucleus is blue by labeling with DAPI. TUNEL assay kit is from Roche and labelled with FITC. Positive apoptosis cells are green.