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Tunel & immunofluorescence double staining protocol (frozen-slides)**1. Apparatus and reagents**

1.1 Apparatus

| Name | Producer | Model |
|---------------------------------|-----------------|---------------------|
| Freezing microtome | Thermo | Cryotome E |
| Glass microscope slides | wanwu | |
| Coverslips | CITOTEST | 10212432C |
| Rocker | wanwu | TSY-B |
| Vortex | wanwu | MX-F |
| Micro-centrifuge | wanwu | D1008E |
| Pipettor | Dragon | KE0003087/KA0056573 |
| Liquid blocker pen | wanwu | WG1066-1 |
| Refrigerator | Haier | BCD-192TGN |
| Ortho-Fluorescent Microscopy | Nikon | NIKON ECLIPSE C1 |
| Imaging system | Nikon | NIKON DS-U3 |

1.2 Major reagents

| Name | Producer | Model |
|--|-----------------|--------------|
| 4% paraformaldehyde | wanwu | G1101 |
| PBS solution | wanwu | G0002 |
| Proteinase K | wanwu | G1205 |
| Permeabilize solution | wanwu | G1204 |
| BSA | wanwu | G5001 |
| Tunel assay kit | | G1501 |
| Primary antibody | | |
| Fluorescent-labelled secondary anyibody | | |
| DAPI | wanwu | G1012 |
| anti-fade mounting medium | wanwu | G1401 |

2.Procedure

2.1 Fix frozen-slides: restore frozen-slides to room temperature and dry in air. Fix in cold acetone 10 min and then dry in air. Wash three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.2 Antigen retrieval: eliminate obvious liquid, mark the objective tissue with liquid blocker pen. Add proteinase K working solution to cover objectives and incubate at 37°C for 25 min. then wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. (Method for configuring working solution of proteinase K, stock solution: PBS=1:9)

2.3 Permeabilization: eliminate excess liquid, add permeabilize working solution to cover objective tissue, then incubate at room temperature for 20 min. wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. (The membrane breaking fluid is 0.1% triton. Configuration method, triton stock solution: PBS=1:1000)

2.4 Equilibrium at room temperature: After the slices are slightly dried, buffer is added to the tissues in the circle, and the buffer is incubated at room temperature for 10 minutes

2.5 Tunel reaction: Take appropriate amount of TDT enzyme, dUTP and buffer in the tunel kit according to the number of slices and tissue size and mix at 1:5:50 ratio. Prepare this reaction solution according to demand before use. Add this mixture to objective tissue placed in a flat wet box, incubate at 37°C for 2 h. be sure to keep the wet box moist by adding water.

2.6 BSA blocking: wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Add 3% BSA to cover the marked tissue to block non-specific binding for 30 min BSA.

2.7 Primary antibody: throw away the blocking solution slightly. Incubate slides with primary antibody (diluted with PBS appropriately) overnight at 4 °C, placed in a wet box containing a little water.

2.8 Secondary antibody: wash slides three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then throw away liquid slightly. Cover objective tissue with secondary antibody (appropriately respond to primary antibody in species), incubate at room temperature for 50 min in dark condition.

2.9 DAPI counterstain in nucleus: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Then incubate with DAPI solution at room temperature for 10 min, kept in dark place.

2.10 Mount: wash three times with PBS (pH 7.4) in a Rocker device, 5 min each. Throw away liquid slightly, then coverslip with anti-fade mounting medium.

2.11 Microscopy detection and collect images by Fluorescent Microscopy. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FITC glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

3.Results

The nuclei stained by DAPI are blue under the excitation of ultraviolet light, the tunel kit is labeled with FITC fluorescein, the positive apoptotic cells are labeled with green light, and the primary antibody is labeled with red.