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<b>Name</b>	<b>Producer</b>	<b>Model</b>
Dehydrator	DIAPATH	Donatello
Glass slide	wanwu	
Coverslip	Citotest	10212432C
Pathologic microtome	Leica	RM2016
Enzyme-free centrifuge tubes	wanwu	EP-150-M
Shaker	wanwu	TSY-B
Vortex	wanwu	MX-F
Pipettor	Dragon	KE0003087/KA0056573
Gene tech pen	Gene tech	GT101
Microscopy	Nikon	NIKON ECLIPSE CI
Imaging system	Nikon	NIKON DS-U3
Incubator	LABOTERY	GSP-70

---

autoclave	PANASONIC	MLS-3751L-PC
Paraffin embedding machine	WHJJ	JB-P5
Frozen flat	WHJJ	JB-L5
Tissue spreader	Kedee	KD-P

---

### 1.2 Major reagents

---

reagent	manufacturer	article number
DEPC	Amresco	E174
4% of paraformaldehyde (DEPC water)	wanwu	G1113
Paraffin wax	Sakura	G1113
Ethanol	SCRC	100092683
Xylene	SCRC	10023418
PBS solution (DEPC)	wanwu	G0020
20×SSC solution	wanwu	G3016-4
BSA	wanwu	G5001
Proteinase K	wanwu	G1205
DAPI	wanwu	G1012
Anti-fluorescence quenching sealing tablets	wanwu	G1401
hybridization buffer	wanwu	G3016-3
anti-DIG-HRP	Jackson	G3016-3
FITC-TSA	wanwu	200-032-156

---

## 2. The steps of the experiment

2.1 **Organization fixed:** take out the organization, wash clean, then Immediately put in the fixed fluid (DEPC) above 12h.

2.2 **Dehydration:** The tissue is dehydrated by gradient alcohol, paraffin, embedding.

2.3 **Section:** the paraffin is sliced through the slicer, the piece of the slice machine and the 62 - degree oven roast for 2 hours.

2.4 **Dewaxing and dehydration:** soak sections in 2 changes of xylene, 15 minutes each. Dehydrate in 2 changes of pure ethanol for 5 minutes each. Then, followed respectively by dehydrating in gradient ethanol of 85%and 75% ethanol 5 minutes for each. Wash in DEPC dilution.

2.5 **Digestion:** according to the tissue fixation time, the slices are boiled in the retrieval solution for 10-15 minutes and naturally cooled. Mark the objective tissue with liquid blocker pen, according to the characteristics of tissues, Add proteinase K(20 ug/ml) working solution to cover objectives and incubate at 37°C for \_\_\_min. Washing with pure water, then wash three times with

PBS (pH 7.4) in a Rocker device, 5 min each.

**2.6 Block endogenous peroxidase:** add 3% methanol-H<sub>2</sub>O<sub>2</sub>, incubate in dark at room temperature for 15min. Wash slides in PBS(PH7.4) three times for 5 min each with gentle agitation.

**2.7 Pre-hybridization:** add Pre-hybridization solution to each section and incubate for 1 h at 37°C.

**2.8 Hybridization:** remove the pre-hybridization solution, add the probe hybridization solution with concentration of \_\_\_\_, and incubate the section in a humidity chamber and hybridize overnight at \_\_\_\_°C.

**2.9 Washing:** remove the hybridization solution. Wash sections with 2×SSC for 10 min at 37°C , Wash sections in 1×SSC two times for 5 min each at 37°C, and wash in 0.5×SSC for 10 min at room temperature. Formamide washing can be added if there are more non-specific hybrids.

**2.10 Blocking:** add blocking solution (rabbit serum) to the section and incubate at room temperature for 30 min.

**2.11** Add the mouse anti-digoxigenin-labeled peroxidase (anti-DIG-HRP): remove the blocking solution and add anti-DIG-HRP. Incubate at 37 °C for 40 min, then Wash sections in PBS four times for 5 min each.

**2.12 TSA developing:** dry sections slightly, add fresh prepared TSA chromogenic reagent to marked tissue. Reaction in dark for 5 min at room temperature. Then wash sections in PBS three times for 10 min each.

**2.13 Stain cell nuclei (counter stain):** incubate with DAPI for 8min in the dark, and then mounting.

**2.14 Microscopic examination and photography :** to take photos with positive fluorescence microscope. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FAM glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

### **3. Interpretation of the results**

The nuclear stained by DAPI were blue under ultraviolet excitation, and the positive expression was a kind of fluorescence labeled by corresponding luciferin. FAM (488) appears green, cy3 appears red. The results of mRNA in situ hybridization were cytoplasmic positive and a few nuclear positive were normal. MicRNA and lncRNA were expressed differently. According to the expression, Different fluorescence brightness is strong or weak.

Note: All reagents, instrument need RNase free.

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HEFEI WANWU technology CO.,LTD