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Cell climbing –DIG -TSA double fluorescent -ISH protocol**1. Apparatus and reagents****1.1 Apparatus**

Name	Producer	Model
Coverslip	Citotest	10212432C
Enzyme-free centrifuge tube	Wanwu	EP-150-M
Shaker	Wanwu	TSY-B
Vortex	Wanwu	MX-F
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Gene tech	GT1001
Refrigerator	HAIER	BCD-192TGN
Microscopy	Nikon	NIKON ECLIPSE CI
Imaging system	Nikon	NIKON DS-U3
Incubator	LABOTERY	GSP-70
Autoclave	PANASONIC	MLS-3751L-PC

1.2 Major reagents

reagent	manufacturer	article number
DEPC	Amresco	E174
4% of paraformaldehyde (DEPC water)	Wanwu	G1113
Ethanol	SCRC	100092683
PBS solution (DEPC)	Wanwu	G0020
20×SSC solution	Wanwu	G3016-4
BSA	Wanwu	G5001
Proteinase K	Wanwu	G1205
DAPI	Wanwu	G1012
Anti-fluorescence quenching sealing tablets	Wanwu	G1401
hybridization buffer	Wanwu	G3016-3
anti-DIG-HRP	Jackson	G3016-3
FITC-TSA	Wanwu	200-032-156

2. The steps of the experiment

2.1. cell climbing fixation: cell climbing was fixed in 4% paraformaldehyde (DEPC) for 20 min, washed 3 times with PBS (pH 7.4) on a decolorizing shaker for 5 min each time.

2.2. Digestion: Mark the objective tissue with liquid blocker pen, according to the characteristics

of different tissues and different indicators, Add proteinase K (20 ug/ml) to cover tissues and incubate at 37°C for _____minutes. washed with sterilized water, then washed three times with PBS, 5 min each time

2.3. **Prehybridization:** add hybridization buffer onto specimen and incubate at 37°C for 1h.

2.4. **Hybridization:** remove the pre-hybridization solution, add the _____ probe hybridization solution with concentration of____, and incubate the cell climbing in a humidity chamber and hybridize overnight at ____°C.

2.5. **Washing:** Remove the hybridization solution. Wash sections with 2×SSC for 5 min at 37°C , wash sections in 1×SSC two times for 5 min each at 37°C, and wash in 0.5×SSC for 10 min at room temperature. Formamide washing can be added if there are more non-specific hybrids.

2.6. **Blocking:** add blocking solution (Rabbit serum) to the section and incubate at room temperature for 30 min.

2.7. Add the mouse anti-digoxigenin-labeled peroxidase (anti-DIG-HRP): Remove the blocking solution and add anti-DIG-HRP. Incubate at 37 °C for 50 min, then Wash sections in PBS three times for 5 min each.

2.8 **FITC-TSA developing:** dry sections slightly, add fresh prepared FITC-TSA chromogenic reagent to marked tissue. Reaction in dark for 5 min at room temperature. Then wash sections in TBST three times for 10 min each and use PBS wash for 5 min.

2.9 **Hybridization:** remove the pre-hybridization solution, add the _____ probe hybridization solution with concentration of____, and incubate the cell climbing in a humidity chamber and hybridize overnight at ____°C.

2.10 **Washing:** Remove the hybridization solution. Wash sections with 2×SSC for 5 min at 37°C , wash sections in 1×SSC two times for 5 min each at 37°C, and wash in 0.5×SSC for 10 min at room temperature. Formamide washing can be added if there are more non-specific hybrids.

2.11 **Blocking:** add blocking solution (Rabbit serum) to the section and incubate at room temperature for 30 min.

2.12 Add the mouse anti-digoxigenin-labeled peroxidase (anti-DIG-HRP): Remove the blocking solution and add anti-DIG-HRP. Incubate at 37 °C for 50 min, then Wash sections in PBS three times for 5 min each.

2.13 **CY3-TSA developing:** dry sections slightly, add fresh prepared CY3-TSA chromogenic reagent to marked tissue. Reaction in dark for 5 min at room temperature. Then wash sections in TBST three times for 10 min each and use PBS wash for 5 min.

2.14 **Stain cell nuclei (counter stain):** Incubate with DAPI for 8min in the dark, and then mounting.

2.15 **Microscopic examination and photography :** To take photos with positive fluorescence microscope. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength

420 nm; FAM glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

3. Interpretation of the results

The nuclear stained by DAPI were blue under ultraviolet excitation, and the positive expression was a kind of fluorescence labeled by corresponding luciferin. FAM (488) appears green, cy3 appears red. The results of mRNA in situ hybridization were cytoplasmic positive and a few nuclear positive were normal. MicRNA and lncRNA were expressed differently. According to the expression, Different fluorescence brightness is strong or weak.

Note: All reagents, instrument need RNase free.

Attached table 1 probe information.