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Frozen Tissue-DIG-TSA-ISH protocol**1. Apparatus and reagents****1.1 Apparatus**

Name	Producer	Model
Frozen slicer	Thermo	Cryotome E
Glass slide	Wanwu	
Coverslip	Citotest	10212432C
Enzyme-free centrifuge tubes	Wanwu	EP-150-M
Shaker	Wanwu	TSY-B
Vortex	Wanwu	MX-F
Pipettor	Dragon	KE0003087/KA0056573
Gene tech pen	Gene tech	GT1001
Refrigerator	HAIER	BCD-192TGN
Microscopy	Nikon	NIKON ECLIPSE CI
Imaging system	Nikon	NIKON DS-U3
Incubator	LABOTERY	GSP-70
autoclave	PANASONIC	MLS-3751L-PC

1.2 Major reagents

reagent	manufacturer	article number
DEPC	Amresco	E174
4% of paraformaldehyde (DEPC water)	Wanwu	G1113
OCT embedding agent	Wanwu	G6059-110ML
Ethanol	SCRC	100092683
Xylene	SCRC	10023418
Sucrose	Wanwu	G5031
PBS solution (DEPC)	Wanwu	G0020
20×SSC solution	Wanwu	G3016-4
BSA	Wanwu	G5001
Proteinase K	Wanwu	G1205
DAPI	Wanwu	G1012
Anti-fluorescence quenching sealing tablets	Wanwu	G1401
anti-DIG-HRP	Jackson	G3016-3
FITC-TSA	Wanwu	200-032-156

2. The steps of the experiment

- 2.1. **Tissue fixation:** tissue was washed and fixed with 4% paraformaldehyde for 12 h or more.
- 2.2. **Tissue dehydration:** Transfer tissue to 15% sucrose solution for 8 h and were subsequently to a 30% sucrose solution overnight.
- 2.3. **Frozen section fixation:** The frozen sections were dried at room temperature, fixed in 4% paraformaldehyde (DEPC) for 10 min, and shaken on the shaker in PBS (pH 7.4) three times for 5 min each.
- 2.4. **Digestion:** Mark the objective tissue with liquid blocker pen, according to the characteristics of different tissues and different indicators, add proteinase K (20 ug/ml) to cover tissues and incubate at 37°C for ___ minutes. Washed with sterilized water and then washed three times with PBS for 5 min each time.
- 2.5. **Block endogenous peroxidase:** add 3% methanol-H₂O₂, incubate in dark at room temperature for 15min. Wash slides in PBS(PH7.4) three times for 5 min each with gentle agitation.
- 2.6. **Prehybridization:** add hybridization buffer onto specimen and incubate at 37°C for 1h.
- 2.7. **Hybridization:** remove the pre-hybridization solution, add the probe hybridization solution with concentration of ___, and incubate the section in a humidity chamber and hybridize overnight at ___°C.
- 2.8. **Washing:** Remove the hybridization solution. Wash sections in 2×SSC for 10 min at 37°C , 1×SSC two times for 5 min each at 37°C, and 0.5×SSC for 10 min at room temperature. Formamide washing can be added if there are more non-specific hybrids.
- 2.9. **Blocking:** add blocking solution (Rabbit serum) to the section and incubate at room temperature for 30 min.
- 2.10. Add the mouse anti-digoxigenin-labeled peroxidase (anti-DIG-HRP): Remove the blocking solution and add anti-DIG-HRP. Incubate at 37 °C for 50 min, then Wash sections in PBS three times for 5 min each.
- 2.11. **TSA developing:** dry sections slightly, add fresh prepared TSA chromogenic reagent to marked tissue. Reaction in dark for 5 min at room temperature. Then wash sections in PBS three times for 5 min each.
- 2.12. **Stain cell nuclei (counter stain):** Incubate with DAPI for 8min in the dark, and then mounting with anti-fluorescence quenching sealing tablets.
- 2.13. **Microscopic examination and photography:** To take photos with positive fluorescence microscope. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FAM glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

3. Interpretation of the results

The nuclear stained by DAPI were blue under ultraviolet excitation, and the positive expression was a kind of fluorescence labeled by corresponding luciferin. FAM (488) appears green, cy3 appears red. The results of mRNA in situ hybridization were cytoplasmic positive and a few nuclear positive were normal. MicRNA and lncRNA were expressed differently. According to the expression, Different fluorescence brightness is strong or weak.

Note: All reagents, instrument need RNase free.

Attached table 1 probe information.