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**1.2 Major reagents**

<b>reagent</b>	<b>manufacturer</b>	<b>article number</b>
Ethanol	SCRC	100092683
Xylene	SCRC	10023418
PBS solution	Wanwu	G0002
Proteinase K	Wanwu	G1205
4% of paraformaldehyde (DEPC water)	Wanwu	G1113
Methanol	SCRC	
Hydrogen peroxide	SCRC	
DEPC	Amresco	E174
20×SSC solution	Wanwu	G3016-4
BSA	Wanwu	G5001
Hematoxylin dye	Wanwu	G1004
neutral balsam	SCRC	
hybridization buffer	Wanwu	G3016-3

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邮箱：2028438226@qq.com

---

anti-DIG-HRP	Jackson	
DAB reagent	Wanwu	G1211

---

## 2. The steps of the experiment

2.1 **Organization fixed:** Take out the organization, Wash clean, then Immediately put in the fixed fluid (DEPC) to fix 2-12h.

2.2 **Dehydration:** The tissue is dehydrated by gradient alcohol, paraffin, Embedding.

2.3 **Section:** The paraffin is sliced through the slicer, the piece of the slice machine and the 62 - degree oven roast for 2 hours.

2.4 **Dewaxing and dehydration:** Soak sections in 2 changes of xylene, 15 minutes each. Dehydrate in 2 changes of pure ethanol for 5 minutes each. Then, followed respectively by dehydrating in gradient ethanol of 85% and 75% ethanol 5 minutes each. Wash in DEPC dilution.

2.5: **Digestion:** According to the tissue fixation time, the slices are boiled in the retrieval solution for 10-15 minutes and naturally cooled. Mark the objective tissue with liquid blocker pen, according to the characteristics of tissues, Add proteinase K(20ug/ml) working solution to cover objectives and incubate at 37°C for \_\_\_min. Washing with pure water, then wash three times with PBS (pH 7.4) in a Rocker device, 5 min each.

2.6 **Block endogenous peroxidase:** add 3% methanol-H<sub>2</sub>O<sub>2</sub> , incubate in dark at room temperature for 15min. Wash slides in PBS(PH7.4) three times for 3 min each, with gentle agitation and protected from light.

2.7 **Pre-hybridization:** Add Pre-hybridization solution to each section and incubate for 1 h at 37°C.

2.8 **Hybridization:** Discard the pre-hybridization solution, add the probe hybridization solution ,concentration\_\_\_, and incubate the section in a humidity chamber and hybridize overnight at \_\_\_°C.

2.9 **Washing:** Remove the hybridization solution. Wash sections with 2×SSC for 5 min at 37°C , Wash sections in 1×SSC two times for 5 min each at 37°C, and wash in 0.5×SSC for 10 min at room temperature. Formamide washing can be added if there are more non-specific hybrids.

2.10 **Blocking:** add blocking solution (Rabbit serum) to the section and incubate at room temperature for 30 min.

2.11 Add the mouse anti-digoxigenin-labeled peroxidase (anti-DIG-HRP): Remove the blocking solution and add anti-DIG-HRP. Incubate at 37 °C for 40 min, then Wash sections in PBS four times for 5 min each.

2.12 **DAB developing:** dry sections slightly, add fresh prepared DAB chromogenic reagent to marked tissue. Manage reaction time by observing in microscopy until positive expression appears brown-yellow. Then stop developing reaction by wash in running tap water.

2.13 Counterstain in nucleus with Hematoxylin staining solution for 3 min and wash in tap water.

---

Treat with the differentiate solution for a few seconds, wash in running tap water. Back to blue by bluing solution, wash in running tap water.

2.14 Dehydrate successively in gradient ethanol of 75%, 85%, and 2 changes of pure ethanol, respectively, 6 min each. Clear in xylene for 6 min and mount with resin mounting medium.

2.15 Image the results using a brightfield fluorescence microscope. Store slides at RT.

### **3. Interpretation of the results**

The nucleus stained with hematoxylin are blue, and the positive expression of DAB was brownish yellow. The results of in situ hybridization of mRNA were cytoplasmic positive and a few nuclear positive were normal. MicRNA and lncRNA were expressed differently. According to the expression of color depth.

Note: All reagents, instrument need RNase free.

**Attached table 1 probe information.**