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**Fluorescence probe-FISH assay protocol (plant paraffin section)****1. Apparatus and reagents****1.1 Apparatus**

<b>Name</b>	<b>Producer</b>	<b>Model</b>
Dehydrator	DIAPATH	Donatello
RNase-free Glass microscope slides	Wanwu	
Coverslip	Citotest	10212432C
Pathologic microtome	Leica	RM2016
Enzyme-free centrifuge tubes	Wanwu	EP-150-M
Shaker	Wanwu	TSY-B
Vortex	Wanwu	MX-F
Pipettor	Dragon	KE0003087/KA0056573
Liquid blocker pen	Gene tech	GT1001
Refrigerator	HAIER	BCD-192TGN
Microscopy	CIC	XSP-C204
Fluorescent Microscopy	NIKON	Nikon Eclipse ci
Imaging system	NIKON	Nikon DS-U3
Incubator	LABOTERY	GSP-70
autoclave	PANASONIC	MLS-3751L-PC
Paraffin embedding machine	WHJJ	JB-P5
Frozen flat	WHJJ	JB-L5
Tissue spreader	Kedee	KD-P

**1.2 Major reagents**

<b>reagent</b>	<b>manufacturer</b>	<b>article number</b>
Ethanol	SCRC	100092683
Xylene	SCRC	10023418
PBS solution	Wanwu	G0002
Proteinase K	Wanwu	G1205
4% of paraformaldehyde (DEPC water)	Wanwu	G1113
DEPC	Amresco	E174
20×SSC solution	Wanwu	G3016-4

---

BSA	Wanwu	G5001
Nuclear fast red dye	Wanwu	G1035
hybridization buffer	Wanwu	G3016-3

---

## 2. Procedure

2.1. **Organization fixed:** Take out the organization, wash clean, then Immediately put in the fixed fluid (DEPC) to fix above 12h.

2.2. **Dehydration:** The tissue is dehydrated by gradient alcohol, paraffin\_ embedding and Vacuum pumping in dehydration process.

2.3. **Section:** The paraffin is sliced through the slicer, the piece of the slice machine and the 62 - degree oven roast for 2 hours.

2.4. **Dewaxing and dehydration:** Soak sections in 2 changes of xylene, 15 minutes each. Dehydrate in 2 changes of pure ethanol for 5 minutes each. Then, followed respectively by dehydrating in gradient ethanol of 85% and 75% ethanol 5 minutes each. Wash in DEPC dilution.

2.5. **Digestion:** Mark the objective tissue with liquid blocker pen, according to the characteristics of tissues, add proteinase K(20 ug/ml) working solution to cover objectives and incubate at 37°C for \_\_\_min. Washing with pure water, then wash three times with PBS (pH 7.4) in a r0ocker device, 5 min each.

2.6. **Pre-hybridization:** Add Pre-hybridization solution to each section and incubate for 1 h at 37°C.

2.7. **Hybridization:** Discard the pre-hybridization solution, add the probe hybridization solution ,concentration\_\_\_\_, and incubate the section in a humidity chamber and hybridize overnight at \_\_\_°C.

2.8. **Washing:** Remove the hybridization solution. Wash sections with 2×SSC for 5 min at 37°C, Wash sections in 1×SSC two times for 5 min each at 37°C, and wash in 0.5×SSC for 10 min at room temperature. Formamide washing can be added if there are more non-specific hybrids.

2.9 **Stain cell nuclei (counter stain):** incubate with DAPI for 8min in the dark, and then mounting.

2.10 **Microscopic examination and photography:** to take photos with positive fluorescence microscope. DAPI glows blue by UV excitation wavelength 330-380 nm and emission wavelength 420 nm; FAM glows green by excitation wavelength 465-495 nm and emission wavelength 515-555 nm; CY3 glows red by excitation wavelength 510-560 nm and emission wavelength 590 nm.

## 3. Interpretation of the results

The positive expression of BCIP/NBT was blue-purple. The results of mRNA in situ hybridization showed that cytoplasm was positive and a few nuclear positive were normal. The expression of micRNA and lncRNA was different. According to the expression, the color is very

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Hefei WANWU technology CO., LTD