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## Scanning Electron Microscopy Report

### 1. Apparatus and Reagents

#### 1.1 Major Apparatus

Name	Producer	Model
Critical Point Dryer	Quorum	K850
Lon Sputtering Apparatus	HITACHI	MC1000
Scanning Electron Microscope	HITACHI	SU8100

#### 1.2 Major Reagents

Name	Producer	Code
Fixative for TEM	Wanwu	G1102
Ethanol	Sinaopharm Group Chemical Reagent Co. LTD	100092183
Isoamyl acetate	Sinaopharm Group Chemical Reagent Co. LTD	10003128
PBS	Wanwu	G0002
OsO <sub>4</sub>	Ted Pella Inc	

### 2. Procedure

**2.1 Harvest tissue block and fixation:** Targeted fresh tissues should be selected to minimize mechanical damage such as pulling, contusion and extrusion. Use a sharp blade to cut and harvest fresh tissue blocks quickly within 1-3 minutes. The area of tissue block should be no more than 3 mm<sup>2</sup>. Wash tissues with PBS gently to remove the blood and hair, etc. Label the target side of tissue (the side you want to observe) by any way, such as making cuts on the opposite side. Make sure to protect tissue blocks, especially the target side, from mechanical damage such as forceps extrusion. The washed tissue blocks are immediately fixed by electron microscopy fixative for 2 hours at room temperature, then transferred into 4°C for preservation and transportation.

**For adherent cell:** Seed cells on a sterile cover glass in a petri dish. Remove the culture medium, then wash slide gently with PBS, followed by adding electron microscopy fixative into petri dish. After fixing for 2 hours at room temperature, transfer the petri dish to 4°C for preservation and transport. Note to avoid sever shock which may result in the cells dropping off from cover glass.

**2.2 Post-fix:** Wash tissue blocks with 0.1 M PB (pH 7.4) for 3 times, 15 min each. Then transfer tissue blocks into 1% OsO<sub>4</sub> in 0.1 M PB (pH 7.4) for 1-2 h at room temperature. After that, wash tissue blocks in 0.1M PB (pH 7.4) for 3 times, 15 min each.

#### 2.3 Dehydrate as followed:

30% ethanol for 15 min;

50% ethanol for 15 min;

70% ethanol for 15 min;

80% ethanol for 15 min;

90% ethanol for 15 min;

95% ethanol for 15 min;

Two changes of 100% ethanol for 15 min;

Finally, isoamyl acetate for 15 min.

**2.4 Drying:** Dry samples with Critical Point Dryer.

**Note: Dry samples and inorganic materials should ignore all the steps above and do the following step directly.**

**2.5 Conductive metal coating:** Specimens are attached to metallic stubs using carbon stickers and sputter-coated with gold for 30s.

**2.6** Observe and take images with scanning electron microscope.